




PRODUCT CATALOG

for pressure-relieving positioning aids




Your specialist for pressure-relieving positioning aids

padcon medical is a company that focuses on the production of medical positioning aids. As a manufactory, with production facilities in Saxony, we have a wide range of products for patient positioning. These are used in hospitals to prevent pressure sores, to treat existing decubitus and to prevent the development of intertrigo or contracture. In addition, we produce custom-made products according to individual specifications to meet the requirements of therapeutic facilities. On the following pages, we present you the custom-made products of our manufactory. Further details can be found on our website at www.padcon-medical.de

At a glance

Table of contents

	03	Production
	04 – 05	Material
01	06 – 08	Small positioning aids
02	09 – 16	Head and neck positioning
03	17 – 22	Leg and heel positioning
04	23 – 26	Arm and hand positioning
05	27 – 29	Lateral positioning
06	30 – 34	Abdominal positioning
07	35 – 37	Back positioning and table rests
08	38 – 41	Fixations and brackets
09	42 – 47	X-ray equipment
10	48 – 51	Other
	52	Custom-made products
	56 – 57	Contact



Regional production according to high quality standards **Made in Germany**

Handmade quality from Saxony

The positioning aids available at padcon medical are manufactured in our manufactory near Leipzig in cooperation with manufacturing companies located in Germany.

The production is done by hand for easy handling: With the help of variable cutting patterns, we create blanks of foam, which are then covered with a hygienic protection or artificial leather cover. The covers are closed with hand-sewn Velcro or zip fasteners. If required, the products are completely hand-sewn and welded. Whether as a roll, half roll, rectangle, cushion or wedge – the shape and function of every positioning cushion meets the medical requirements of patient care.

In addition to arm and leg positionings, our wide range of products also includes lateral, abdominal and back positioning sets that meet all normative and regulatory requirements for Class I medical devices according to quality management and carry the CE mark.

All materials used have been tested by accredited test institutes for cytotoxicity and biocompatibility and can therefore be used on patients without restrictions.

The processed covers are not only free of harmful substances and allergens, they can also be reprocessed by means of a standardised wipe disinfection.

It's the correct care that matters

Instructions

Positioning aids with hygienic protection cover

The positioning aids produced at padcon medical are covered with our PU cover in navy blue as standard. This is a hyper-elastic, slip-resistant hygienic protection cover.

In order not to jeopardise the durability of the material, please follow the special washing and care instructions. Washing and care instructions Please use only DGHM-approved disinfectants for normal cleaning and disinfection of PU covers. Polyurethane protection covers can be spray disinfected or wiped with a cloth soaked in disinfectant. The use of disinfectants containing alcohol should be avoided. Coloured disinfectants, e.g. those containing iodine, can cause discolouration, while chemical cleaning agents, e.g. those containing chlorine, cause brittle hardening, which leads to leakage.

Hygienic protection covers such as the PU cover can be washed at 60° according to RAL-GZ 992-3. In order to avoid fluffing of detached fibres and floating particles, we recommend closing all closures before starting the washing programme. We also advise against exposing the covers of the positioning aids to increased heat.

Please note:

Hygienic protection covers may be washed at up to 60°C.

Hygienic protection covers must not be dried in a tumble dryer.

Hygienic protection covers must not be sterilised in the VDV chamber or autoclave.

Hygienic protection covers must not be pressed.

Hygienic protection covers must not be ironed.

As an option, we offer to produce the positioning aids with a Hypro-Elastic cover, washable at 95°, made of slip-resistant polyurethane or artificial leather.

Please contact us for an individual offer.



Standard cover navy blue

PU cover (polyurethane, weldable, slip-resistant)




Special colours


Hypro-Elastic (polyurethane, washable at up to 95°C)



N° 15 – Orange



N° 16 – Green



N° 17 – Creme



N° 18 – Blue



SMALL
POSITIONING AIDS

01

RECTANGLE

with velcro fastener

Dimensions	Item number
400 x 240 x 50 mm	PD00189
400 x 200 x 50 mm	PD00190
300 x 250 x 55 mm	PD00192
300 x 250 x 65 mm	PD00193
300 x 200 x 55 mm	PD00194
300 x 200 x 65 mm	PD00195
240 x 200 x 50 mm	PD00197

Cover · PU slip-resistant navy blue

Core · Spacer material/foam RG 30



Special colours

Cover · Hypro-Elastic (polyurethane)



ROLL

with velcro fastener or with zip fastener

Dimensions	Item number
400 x Ø 150 mm	PD00220
200 x Ø 60 mm	PD00221
200 x Ø 150 mm	PD00222
200 x Ø 180 mm	PD00223
450 x Ø 80 mm	PD00224
500 x Ø 100 mm	PD00225
300 x Ø 100 mm	PD00226
300 x Ø 150 mm	PD00227
120 x Ø 55 mm	PD00228
500 x Ø 120 mm	PD01031*
500 x Ø 150 mm	PD01034*

Cover · PU slip-resistant navy blue

Core · Spacer material/viscoelastic foam/*foam RG 30



Special colours

Cover · Hypro-Elastic (polyurethane)



BABY POSITIONING ROLL

completely sewn

Dimensions	Item number
230 x Ø 40 mm	PD07710



Special colours

Cover · Hypro-Elastic (polyurethane)



Cover · PU slip-resistant navy blue

Core · Microbeads

HALF ROLL

with velcro fastener

Dimensions	Item number
220 x 230 x 115 mm	PD00229
400 x 230 x 115 mm	PD00230
500 x 230 x 115 mm	PD00231
600 x 230 x 115 mm	PD00232
500 x 200 x 100 mm	PD00233

Cover · PU slip-resistant navy blue

Core · Spacer material/viscoelastic foam



Special colours

Cover · Hypro-Elastic (polyurethane)



THREE-QUARTER ROLL

with velcro fastener or with zip fastener on three sides

Dimensions	Item number
250 x 120 x 90 mm	PD00240
500 x 220 x 165 mm	PD00241*
300 x 140 x 110 mm	PD00242*
150 x 120 x 90 mm	PD00243
500 x 150 x 100 mm	PD00245*
500 x 100 x 75 mm	PD00246*
400 x 150 x 110 mm	PD01010

Cover · PU slip-resistant navy blue

Core · Spacer material/viscoelastic foam/*foam RG 30



Special colours

Cover · Hypro-Elastic (polyurethane)



WEDGE

with velcro fastener or with zip fastener

Dimensions	Item number
400 x 400 x 240/0 mm Wedge 30°	PD01440
400 x 600 x 200/0 mm Wedge 18°	PD01450
400 x 800 x 220/0 mm Wedge 15°	PD01460

Cover · PU slip-resistant navy blue

Core · foam RG 30



Special colours

Cover · Hypro-Elastic (polyurethane)





HEAD AND NECK
POSITIONING

02

HEAD POSITIONING CUSHION

with velcro fastener

Dimensions	Item number
250 x 300 x 100/60 mm	PD00270
500 x 300 x 100/60 mm	PD00271
250 x 300 x 100/60/70 mm	PD00273
500 x 300 x 100/60/70 mm	PD00274



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/viscoelastic foam

ITS HEAD POSITIONING CUSHION

with velcro fastener

Dimensions	Item number
340 x 250 x 100/55/100 mm	PD00280



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/viscoelastic foam

HEAD POSITIONING CUSHION

completely sewn

Dimensions	Item number
250 x 250 x 110/60 mm	PD00450
250 x 250 x 55/30 mm	PD00451



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/foam RG 30

HEAD RING

with velcro fastener

Dimensions	Item number
O Ø 230, I Ø 100 Height 50 mm	PD00260



Cover • PU slip-resistant navy blue
Core • Spacer material/viscoelastic foam

Special colours    
Cover • Hypro-Elastic (polyurethane)



HEAD RING ROUND PLATE

completely sewn/round plate O Ø 250 mm, Height 15 mm

Dimensions	Item number
O Ø 260, I Ø 100 Height 50 mm	PD00265
O Ø 260, I Ø 100 Height 70 mm	PD00266



Cover • PU slip-resistant navy blue
Core • Spacer material/foam RG 30

Special colours    
Cover • Hypro-Elastic (polyurethane)



HEAD RING WITHOUT ROUND PLATE

completely sewn

Dimensions	Item number
O Ø 260, I Ø 100 Height 50 mm	PD00267



Cover • PU slip-resistant navy blue
Core • Spacer material/foam RG 30

Special colours    
Cover • Hypro-Elastic (polyurethane)



BABY HEAD RING

completely sewn

Dimensions	Item number
Ø 120, I 50 mm	PD07720
Ø 150, I 70 mm	PD07721



Cover • PU slip-resistant navy blue
Core • Microbeads

Special colours    

Cover • Hypro-Elastic (polyurethane)

TODDLER HEAD RING

completely sewn

Dimensions	Item number
Ø 170, I 80 mm	PD07722



Cover • PU slip-resistant navy blue
Core • Microbeads

Special colours    

Cover • Hypro-Elastic (polyurethane)

BABY POSITIONING CRESTS

completely sewn

Dimensions	Item number
Ø 100, I 30 mm	PD07730
Ø 140, I 40 mm	PD07731



Cover • PU slip-resistant navy blue
Core • Microbeads

Special colours    

Cover • Hypro-Elastic (polyurethane)

HEAD LATERAL POSITIONING CUSHION, left

with zip fastener

Dimensions	Item number
140 x Ø 300 mm	PD01500
100 x Ø 300 mm flat	PD01501



Cover • PU slip-resistant navy blue
Core • Spacer material/foam RG 35

Special colours    
Cover • Hypro-Elastic (polyurethane)



HEAD CUSHION

with velcro fastener or with zip fastener

Dimensions	Item number
400 x 400 x 100 mm	PD00282
200 x 200 x 50 mm	PD00283
200 x 300 x 40 mm	PD01050



Cover • PU slip-resistant navy blue
Core • foam RG 20

Special colours    
Cover • Hypro-Elastic (polyurethane)



NOSE CUSHION

with velcro fastener

Dimensions	Item number
320 x 320 x 70 mm	PD00200
320 x 320 x 140 mm	PD00210



Cover • PU slip-resistant navy blue
Core • foam RG 35/viscoelastic foam

Special colours    
Cover • Hypro-Elastic (polyurethane)



FACE CUSHION

with zip fastener

Dimensions		Item number
140 x Ø 300 mm	2 indentation top	PD00500
140 x Ø 300 mm	2 indentation top, 2 indentation bottom	PD00501



Cover • PU slip-resistant navy blue
Core • Spacer material/foam RG 35

Special colours    
Cover • Hypro-Elastic (polyurethane)



NECK CUSHION

with velcro fastener

Dimensions	Item number
250 x 300 x 70 / 30 mm	PD00275
500 x 300 x 70 / 30 mm	PD00276
400 x 300 x 70 / 30 mm	PD00277



Cover • PU slip-resistant navy blue
Core • Spacer material/viscoelastic foam

Special colours    
Cover • Hypro-Elastic (polyurethane)



NECK AND ABDOMINAL POSITIONING CUSHION, elliptical

with velcro fastener

Dimensions	Item number
350 x 180 x 90 mm	PD00198



Cover • PU slip-resistant navy blue
Core • Spacer material/viscoelastic foam


Special colours    
Cover • Hypro-Elastic (polyurethane)



HEAD AND ARM POSITIONING CUSHION, oval

with velcro fastener

Dimensions	Item number
350 x 180 x 90 mm	PD00199



Cover · PU slip-resistant navy blue
Core · Spacer material/viscoelastic foam

Special colours


Cover · Hypro-Elastic (polyurethane)

CE 150 60 X X X X X

HEAD AND ARM POSITIONING CUSHION

small with velcro fastener / big with zip fastener on 3 sides

Dimensions	Item number
300 x 250 x 100/50 mm	PD00250
300 x 500 x 100/50 mm	PD00251



Cover · PU slip-resistant navy blue
Core · Spacer material/viscoelastic foam

Special colours

Cover · Hypro-Elastic (polyurethane)

CE 150 60 X X X X X

HEAD AND ARM POSITIONING WEDGE

with velcro fastener

Dimensions	Item number
250 x 300 x 100/50 mm	PD00255



Cover · PU slip-resistant navy blue
Core · Spacer material/viscoelastic foam

Special colours

Cover · Hypro-Elastic (polyurethane)

CE 150 60 X X X X X

SHOULDER REST CUSHION

with zip fastener

Dimensions	Item number
500 x 360 x 160 mm	PD00470



Cover • PU slip-resistant navy blue
Core • Spacer material/ foam RG 30/ Neopolen

Special colours    


Cover • Hypro-Elastic (polyurethane)

INTUBATION CUSHION

with velcro fastener

Dimensions	Item number
300 x 250 x 180/130 mm	PD00370
300 x 250 x 105/60 mm	PD00371



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours    

Cover • Hypro-Elastic (polyurethane)



LEG AND HEEL
POSITIONING

03

HEEL CUSHION

with velcro fastener and indentation

Dimensions	Item number
200 x 70 x 105 mm	PD00288



Cover • PU slip-resistant navy blue
Core • foam RG 56

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

60

HEEL CUSHION

with velcro fastener

Dimensions	Item number
250 x 140 x 80 / 30 / 80 mm	PD00289



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

60

HEEL CUSHION

with velcro fastener

Dimensions	Item number
250 x 200 x 100 / 40 / 100 mm	PD00290



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

60

HEEL CUSHION

with velcro fastener

Dimensions	Item number
500 x 200 x 100 / 40 / 100 mm	PD00292



Cover • PU slip-resistant navy blue
Core • Spacer material / viscoelastic foam

Special colours    

Cover • Hypro-Elastic (polyurethane)

HEEL CUSHION

with velcro fastener and middle raise

Dimensions	Item number
500 x 200 x 100 / 40 / 100 / 40 / 100 mm	PD00293



Cover • PU slip-resistant navy blue
Core • Spacer material / viscoelastic foam

Special colours    

Cover • Hypro-Elastic (polyurethane)

HEEL CUSHION WITH FIXATION

with velcro fastener and velcro fixation

Dimensions	Item number
250 x 200 x 100 / 40 / 100 mm	PD00291



Cover • PU slip-resistant navy blue
Core • Spacer material / viscoelastic foam

Special colours    

Cover • Hypro-Elastic (polyurethane)

NECK AND HEEL CUSHION

with velcro fastener

Dimensions	Item number
250 x 300 x 70 / 30 mm	PD00275
500 x 300 x 70 / 30 mm	PD00276
400 x 300 x 70 / 30 mm	PD00277



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material / viscoelastic foam

LEG CUSHION

with velcro fastener

Dimensions	Item number
250 x 400 x 100 / 40 / 100 mm	PD00285



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material / viscoelastic foam

PELVIC SUPPORT CUSHION

with velcro fastener and velcro / fleece fixation

Dimensions	Item number
230 x 130 x 50 mm	PD00300



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material / foam RG 30

KNEE AND FOOT REST

with velcro fastener

Dimensions	Item number
500 x 550 x 70 / 30 mm	PD00375
500 x 750 x 140 / 60 mm	PD00376



Cover • PU slip-resistant navy blue
Core • Spacer material / viscoelastic foam


Special colours    
 Cover • Hypro-Elastic (polyurethane)



REST FOR GRÖBELS' LEG SUPPORT

with velcro fastener and pocket and velcro fixation

Dimensions	Item number
300 x 450 x 15 mm	PD00180
300 x 550 x 15 mm	PD00181



Cover • PU slip-resistant navy blue
Core • Spacer material / viscoelastic foam

Special colours    
 Cover • Hypro-Elastic (polyurethane)



LEG ELEVATION POSITIONING WEDGE 45°

with zip fastener

Dimensions	Item number
600 x 200 x 220 mm	PD01100
720 x 200 x 220 mm	PD01101



Cover • PU slip-resistant navy blue
Core • foam RG 30


Special colours    
 Cover • Hypro-Elastic (polyurethane)



LEG ELEVATION POSITIONING WEDGE 45°, double

with zip fastener

Dimensions	Item number
600 x 420 x 200 mm	PD01110
750 x 420 x 200 mm	PD01111



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    
Cover • Hypro-Elastic (polyurethane)



LEG POSITIONING WEDGE WITH RECESS

with zip fastener

Dimensions	Item number
800 x 220 x 260 / 50 mm	PD01200
600 x 220 x 200 / 50 mm	PD01205



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    
Cover • Hypro-Elastic (polyurethane)



SPREADING WEDGE WITH RECESS

with zip fastener

Dimensions	Item number
380 x 300 / 150 x 200 mm	PD01120
450 x 360 / 100 x 160 mm	PD01122



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    
Cover • Hypro-Elastic (polyurethane)






ARM REST

with velcro fastener

Dimensions	Item number
------------	-------------

150 x 450 x 15 mm	
without velcro fixation	PD00170
with velcro fixation	PD00171
with velcro fixation and flap	PD00172
with two velcro fixation	PD00173
with two velcro fixation and flap	PD00174



Special colours

Cover • Hypro-Elastic (polyurethane)



ARM REST

with velcro fastener

Dimensions	Item number
------------	-------------

170 x 500 x 15 mm	
with velcro fixation and flap	PD00175



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue


Core • Spacer material/viscoelastic foam

HEAD AND ARM POSITIONING CUSHION

small with velcro fastener / big with zip fastener on 3 sides

Dimensions	Item number
------------	-------------

300 x 250 x 100/50 mm	PD00250
300 x 500 x 100/50 mm	PD00251



Special colours

Cover • Hypro-Elastic (polyurethane)




Cover • PU slip-resistant navy blue

Core • Spacer material/viscoelastic foam

HEAD AND ARM POSITIONING WEDGE

with velcro fastener

Dimensions	Item number
250 x 300 x 100/50 mm	PD00255



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)



ARM POSITIONING CUSHION WITHOUT RECESS

with zip fastener

Dimensions	Item number
200 x 450 x 270/50 mm	PD01201



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)



ARM POSITIONING CUSHION WITH RECESS

with zip fastener

Dimensions	Item number
220 x 650 x 270/30 mm	PD01202
200 x 450 x 270/50 mm	PD01203



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours


Cover • Hypro-Elastic (polyurethane)



INJECTION HAND ROLL

completely sewn

Dimensions	Item number
120 x Ø 55 mm	PD00228



Cover • PU slip-resistant navy blue
Core • Spacer material/ ethafoam

Special colours    

Cover • Hypro-Elastic (polyurethane)

INJECTION CUSHION WITH RECESS

with zip fastener

Dimensions	Item number
150 x 300 x 100 / 30 mm	PD01210



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    

Cover • Hypro-Elastic (polyurethane)



LATERAL POSITIONING

05

LATERAL POSITIONING CUSHION

with zip fastener

Dimensions	Item number
500 x 720 x 150 / 30 mm	PD00155



Cover • PU slip-resistant navy blue

Core • Spacer material / foam RG 30 / viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)



LATERAL POSITIONING CUSHION, rounded

with velcro fastener

Dimensions	Item number
500 x 250 x 100 / 50 / 0 mm	PD00156



Cover • PU slip-resistant navy blue

Core • Spacer material / foam RG 35

Special colours


Cover • Hypro-Elastic (polyurethane)



ITS LATERAL POSITIONING WEDGE

with zip fastener

Dimensions	Item number
600 x 300 x 150 mm	PD00159



Cover • PU slip-resistant navy blue

Core • Spacer material / foam RG 35 / viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)



ITS LATERAL POSITIONING SET 135°, two-part

with velcro fastener / zip fastener

Dimensions	Item number
ITS lateral positioning wedge 600 x 300 x 150 mm	PD00158
Nose cushion, small 320 x 320 x 70 mm	PD00200



Cover • PU slip-resistant navy blue

Core • Spacer material / foam RG 35 / viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)



LATERAL POSITIONING WEDGES 30°, left

with zip fastener

Dimensions	Item number
600 x 300 / 400 x 120 mm	PD01480



Cover • PU slip-resistant navy blue

Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)



LATERAL POSITIONING WEDGES 30°, right

with zip fastener

Dimensions	Item number
600 x 300 / 400 x 120 mm	PD01481



Cover • PU slip-resistant navy blue

Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)





For the treatment of acute gas exchange disorders such as ARDS (Acute Respiratory Distress Syndrome)

ITS abdominal positioning set (PD00149RV)



Application

With the help of a spiral sheet, the patient is turned to the side so that the positioning cushions – the upper and lower part – can be placed underneath. Positioning should always be carried out with three to four nurses. The patient's head is positioned on the face cushion. Thus, eyes, nose and chin are exposed and the respiratory tract is free. Depending on the patient, the arms are positioned alternately upwards and/or downwards for the duration of the positioning. If necessary, foam blanks or absorbent compresses can be placed underneath vulnerable areas of the face to additionally reduce the contact pressure and prevent sticking. If necessary, the thorax can be additionally padded with hydrocolloidal dressing material.

Positioning duration*

16 – 20 hours

Description

ITS abdominal positioning set
special high N, welded seams

Core

Spacer material, foam RG 35,
viscoelastic foam

Cover

PU cover navy blue
with zip fastener

Dimensions

Head part Ø 350x240 mm

Upper part 350x600x 240 mm

Lower part 500 x 600 x 240 mm

*The individual clinical picture must be taken into account!

OPERATION ABDOMINAL POSITIONING SET N, three-part

seams welded/sewn, with zip fastener

Dimensions

Item number

Upper part
310 x 500 x 120 mm
Lower part
310 x 500 x 120 mm
Knee positioning plate
250 x 500 x 25 mm

PD00141RV



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/ foam RG 35

OPERATION ABDOMINAL POSITIONING SET W, three-part

seams welded/sewn, with zip fastener

Dimensions

Item number

Upper part
300 x 500 x 140 mm
Lower part
300 x 500 x 140 mm
Knee positioning plate
250 x 500 x 25 mm

PD00143RV



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/ foam RG 35/
viscoelastic foam

ITS ABDOMINAL POSITIONING SET N, three-part

seams welded/sewn, with zip fastener

Dimensions

Item number

Head part Ø 350 x 240 mm
Upper part 350 x 600 x 240 mm
Lower part 500 x 600 x 240 mm

PD00149RV



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • Spacer material/ foam RG 35/
viscoelastic foam

THORAX POSITIONING CUSHION

with velcro fastener

Dimensions	Item number
600 x 400 x 100 mm	PD00160



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours    

Cover • Hypro-Elastic (polyurethane)

THORAX POSITIONING WEDGE

with zip fastener

Dimensions	Item number
500 x 600 x 100/0 mm	PD00164



Cover • PU slip-resistant navy blue
Core • Spacer material/ foam RG 30

Special colours    

Cover • Hypro-Elastic (polyurethane)

THORAX POSITIONING CUSHION, two-part

with velcro fastener

Dimensions	Item number
770/630 x 210/80 x 50 mm	PD00165
770/630 x 210/80 x 90 mm	PD00166



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours    


Cover • Hypro-Elastic (polyurethane)

POSITIONING CUSHION WITH RECESS

completely sewn

Dimensions	Item number
410 x 720 x 100 mm	PD00168
460 x 820 x 100 mm	PD00169



Cover • PU slip-resistant navy blue
Core • foam RG 56

Special colours    

Cover • Hypro-Elastic (polyurethane)

BREAST CUSHION

with zip fastener and cover flap

Dimensions	Item number
500 x 220 x 245 mm	PD00167



Cover • PU slip-resistant navy blue
Core • foam RG 30 / polystyrene

Special colours    

Cover • Hypro-Elastic (polyurethane)



BACK POSITIONING
AND TABLE RESTS

07

TABLE TOP FOR OPERATING TABLE

with velcro fastener

Dimensions	Item number
960 x 500 x 25 mm	PD00100
1900 x 500 x 25 mm	PD00101
1200 x 500 x 25 mm	PD00106

Cover · Hypro-Elastic (polyurethane)
Underside optionally with slip-resistant
PU as anti-slip cover

Core · Spacer material/ viscoelastic foam



Colours



Cover · Hypro-Elastic (polyurethane)



TABLE TOP SET, three-part

with velcro fastener and velcro fixation

Dimensions	Item number
1 Stück: 960 x 500 x 25 mm,	PD00110
2 Stück: 840 x 240 x 25 mm	

Cover · Hypro-Elastic (polyurethane)
Underside optionally with slip-resistant
PU as anti-slip cover

Core · Spacer material/ viscoelastic foam



Colours



Cover · Hypro-Elastic (polyurethane)



TABLE TOP FOR EXTENSION TABLE

with velcro fastener

Dimensions	Item number
1000 x 500 x 25 mm	PD00102
1150 x 500 x 25 mm	PD00103

Cover · Hypro-Elastic (polyurethane)
Underside optionally with slip-resistant
PU as anti-slip cover

Core · Spacer material/ viscoelastic foam



Colours



Cover · Hypro-Elastic (polyurethane)



TABLE TOP FOR GYNAECOLOGY

with velcro fastener

Dimensions	Item number
1080 x 500 x 25 mm	PD00130



Cover • Hypro-Elastic (polyurethane)
Underside optionally with slip-resistant
PU as anti-slip cover
Core • Spacer material/ viscoelastic foam

Colours

Cover : Hypro-Elastic (polyurethane)

thane)

CE


150

95

BACK SUPPORT CUSHION

with velcro fastener

Dimensions	Item number
500 x 400 x 10/50/10 mm	PD00330



Cover • PU slip-resistant navy blue
Core • Spacer material/ viscoelastic foam

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

60



dreapunk

ARM CUFF

with velcro fastener

Dimensions

250 x 520 x 15 mm

Item number

PD00310



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • PU slip-resistant navy blue

Core • foam RG 20

ARM SUPPORT CUFF FOR ANAESTHESIA BRACKET

padded / velcro / fleece fastener

Dimensions

450 x 180 mm

Item number

PD00311



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • Strong foil / PU slip-resistant navy blue

ARM SUPPORT CUFF

velcro / fleece fastener

Dimensions

400 x 120 mm

Item number

PD00315



Special colours

Cover • Hypro-Elastic (polyurethane)



Cover • Strong foil / PU slip-resistant navy blue

Core • foam RG 30

SUPPORT BELT

velcro / fleece fastener

Dimensions	Item number
575 x 60 mm	PD00320
860 x 55 mm	PD00321
680 x 60 mm	PD00323



Special colours



Cover • Hypro-Elastic (polyurethane)

Cover • Strong foil / PU slip-resistant navy blue



LEG SUPPORT BELT

velcro / fleece fastener

Dimensions	Item number
1750 x 120 x 5 mm	PD00340
1900 x 120 x 5 mm	PD00341



Special colours



Cover • Hypro-Elastic (polyurethane)

Cover • Strong foil / PU slip-resistant navy blue
Core • foam RG 30



KNEE FIXATION CUFF

with two fixations for left and right

Dimensions	Item number
450 x 650 / 470 mm	PD00345



Special colours



Cover • Hypro-Elastic (polyurethane)

Cover • Strong foil / PU slip-resistant navy blue
Core • viscoelastic foam



LEATHER FIXATION BELT

padded / 50 mm round buckle / waxed bridle leather

Dimensions	Item number
1900 x 80 x 4 mm	PD00350
2200 x 120 x 4 mm	PD00351



Cover · Upholstery sleeve, imitation leather
Material · Waxed natural leather



FIXATION BELT, two-part

velcro / fleece fastener

Dimensions	Item number
650 x 120 mm	PD00360 65 N
850 x 120 mm	PD00361 85 N



Cover · PU slip-resistant navy blue
Core · Spacer material

Special colours
 Cover · Hypro-Elastic (polyurethane)



SUPPORT BELT WITH RING

padded / with retaining rings


Dimensions	Item number
520 x 120 mm	PD00365



Cover · Strong foil / PU slip-resistant navy blue
Core · Spacer material / foam RG 30

Special colours
 Cover · Hypro-Elastic (polyurethane)






dropout

X-RAY SET, ten-part

completely welded / sewn

Dimensions		Item number
1 round plate	(PD06001)	PD06000
1 ring plate	(PD06003)	
1 roll	(PD06004)	
1 rectangle	(PD06006)	
2 triangles	(PD06010)	
2 triangles	(PD06012)	
2 wedges	(PD06016)	

Cover · PU slip-resistant navy blue
Core · foam RG 30



Special colours    
Cover · Hypro-Elastic (polyurethane)

CE      

X-RAY SET, nine-part

completely welded / sewn

Dimensions		Item number
1 big triangle	(PD06021)	PD06020
2 medium triangles	(PD06022)	
2 small triangles	(PD06023)	
2 wedges	(PD06024)	
1 rectangle	(PD06025)	
1 rectangle	(PD06026)	

Cover · PU slip-resistant navy blue
Core · foam RG 30



Special colours    
Cover · Hypro-Elastic (polyurethane)

CE      

ROUND PLATE

completely welded / sewn

Dimensions	Item number
Ø 230 x 50 mm	PD06001



Cover · PU slip-resistant navy blue
Core · foam RG 30

Special colours    
Cover · Hypro-Elastic (polyurethane)

CE      

RING PLATE

completely welded / sewn

Dimensions	Item number
Ø 230 x 50 mm	PD06003



Cover · PU slip-resistant navy blue
Core · foam RG 30

Special colours    
Cover · Hypro-Elastic (polyurethane)



ROLL

completely welded / sewn

Dimensions	Item number
Ø 100 x 250 mm	PD06004
Ø 120 x 500 mm	PD06005



Cover · PU slip-resistant navy blue
Core · foam RG 30

Special colours    
Cover · Hypro-Elastic (polyurethane)



RECTANGLE

completely welded / sewn

Dimensions	Item number
180 x 240 x 50 mm	PD06006
360 x 240 x 30 mm	PD06007
360 x 240 x 50 mm	PD06008



Cover · PU slip-resistant navy blue
Core · foam RG 30

Special colours    
Cover · Hypro-Elastic (polyurethane)



TRIANGLE 45°, isosceles

completely welded / sewn

Dimensions	Item number
180 x 240 x 120 / 0 mm	PD06010
180 x 100 x 50 / 0 mm	PD06012



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

WEDGE

completely welded / sewn

Dimensions		Item number
180 x 240 x 70 / 10 mm	Wedge 20°	PD06016
250 x 200 x 70 / 0 mm	Wedge 15°	PD06017
230 x 180 x 100 / 0 mm	Wedge 30°	PD06018
178 x 280 x 76 / 0 mm	Wedge 15°	PD06019
180 x 240 x 60 / 10 mm	Wedge	PD06024



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)


CE

150

TRIANGLE

completely welded / sewn

Dimensions	Item number
big triangle	PD06021
360 x 240 x 110 / 20 mm	
medium triangle	PD06022
180 x 240 x 110 / 20 mm	
small triangle 45 °, isosceles	PD06023
180 x 100 x 50 / 0 mm	



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

150

X-RAY TABLE TOP

completely welded / sewn

Dimensions	Item number
2000 x 600 x 20 mm	PD06142
2000 x 600 x 30 mm	PD06143



Cover • Hypro-Elastic
with optional PU slip-resistant navy blue
Core • foam RG 40

Colours    

Cover • Hypro-Elastic (polyurethane)

TABLE TOP FOR MRT

completely welded / sewn

Dimensions	Item number
1500 x 400 x 10 / 40 / 10 mm	PD06200



Cover • Hypro-Elastic
with optional PU slip-resistant navy blue
Core • foam RG 30

Colours    

Cover • Hypro-Elastic (polyurethane)

X-RAY HEAD BOWL

completely welded / sewn

Dimensions	Item number
250 x 250 x 100 mm	PD06330



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    

Cover • Hypro-Elastic (polyurethane)

LEG RECESS


completely welded / sewn

Dimensions

150 x 180 x 80 mm
180 x 210 x 80 mm

Item number

PD06350
PD06351



Cover • PU slip-resistant navy blue
Core • foam RG 56

Special colours

Cover • Hypro-Elastic (polyurethane)



LEG POSITIONING CUSHION FOR MRT, two coils recess


completely welded / sewn

Dimensions

550 x 400 x 170 mm

Item number

PD06360



Cover • PU slip-resistant navy blue
Core • foam RG 56

Special colours

Cover • Hypro-Elastic (polyurethane)



LEG POSITIONING CUSHION FOR MRT, underside convexed

completely welded / sewn

Dimensions

440 x 730 x 230 mm

Item number

PD06361



Cover • PU slip-resistant navy blue
Core • foam RG 40

Special colours

Cover • Hypro-Elastic (polyurethane)





TUNNEL CUSHION

with zip fastener

Dimensions	Item number
250 x 300 x 280 mm	PD00530
350 x 500 x 220 mm	PD00531



Cover • PU slip-resistant navy blue
Core • Spacer material/ foam RG 30/ ethafoam

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

BONE CUSHION

with zip fastener

Dimensions	Item number
300 x 140 x 140 mm	PD01305
260 x 110 x 110 mm	PD01306



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

VEIN CUSHION

with zip fastener

Dimensions	Item number
600 x 450 x 100 mm	PD01330
550 x 750 x 200 mm	PD01331



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

SEAT CUSHION FOR CARE CHAIR

with zip fastener

Dimensions	Item number
550 x 490 x 60 mm	PD01434



Cover • PU slip-resistant navy blue
Core • foam RG 56 / Spacer material

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

SEAT CUSHION, round

with zip fastener

Dimensions	Item number
O Ø 400, I Ø 190, Height 60 mm	PD01435



Cover • PU slip-resistant navy blue
Core • Spacer material / foam RG 56

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

SEAT CUSHION, round, front-open

with zip fastener

Dimensions	Item number
O Ø 480, I Ø 170, Height 90 mm	PD01438



Cover • PU slip-resistant navy blue
Core • foam HR 43

Special colours

Cover • Hypro-Elastic (polyurethane)

CE

LOWER EXTREMITY POSITIONING CUSHION

with zip fastener

Dimensions	Item number
600 x 400 x 300 mm	PD01340
500 x 600 x 400 mm	PD01341
400 x 500 x 300 mm	PD01342



Cover • PU slip-resistant navy blue
Core • foam RG 30

Special colours    

Cover • Hypro-Elastic (polyurethane)

CUBE

with zip fastener

Dimensions	Item number
300 x 300 x 300 mm	PD01370
400 x 400 x 400 mm	PD01371



Cover • PU slip-resistant navy blue
Core • Fester Core

Special colours    

Cover • Hypro-Elastic (polyurethane)

Individual design possibilities

Custom-made products

Since we produce positioning aids by hand, we are open to individual wishes. We are happy to produce custom-made products according to your needs.

Abdominal positioning cushion for spinal operations



Heel free positioning ring



Arm support cuff (double)



Rescue head positioning cushion



Knee cuff padded (for obese patients)



We produce according to high
quality standards in Germany
- **but in the name of health,
we deliver without borders.**

Get to know padcon medical

the manufacturer for positioning aids

+49 (0) 34299 - 704160
info@padcon-medical.de

You can reach us Monday to Friday from 9:00 to 16:00 o'clock.

Do you have any questions for padcon medical?

We would be pleased to advise you in detail about our product variety and the associated production processes. From preparation to custom-made products, we implement your orders reliably and solution-oriented. In addition, as experts in positioning aids, we offer a wide range of training courses for the correct use of our positioning cushions. Contact us.

Further details can be found on our website.
www.padcon-medical.de

Imprint

Publisher padcon medical GmbH
Hauptstraße 16
04564 Böhlen -DE
T +49 (0) 34299 - 704160
F +49 (0) 34299 - 704161
info@padcon-medical.de
www.padcon-medical.de

Concept dreipunkt
Text Wallberg & Zúñiga GbR
Layout Körnerstrasse 56
Photography 04107 Leipzig - DE
T +49 (0) 341 - 14990310
F +49 (0) 341 - 14990311
info@dreipunkt.design
www.dreipunkt.design

Image credits padcon medical GmbH
3, 7, 8, 10-16, 18-22, 24-26, 28, 29, 32-34, 36, 37, 39-41, 43-47, 49-51, 52, 54, 55
dreipunkt
Cover, 1, 4, 5, 6, 7, 8, 9, 14, 15, 17, 18, 19, 20, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 35, 36, 37, 38,
40, 42, 46, 48, 50, 51, 53



padcon medical GmbH

Hauptstraße 16
04564 Böhlen – DE

T +49 (0) 34299 – 704160
F +49 (0) 34299 – 704161

info@padcon-medical.de
www.padcon-medical.de